

## Moisture Analyser BM-60

- new, state of the art moisture analyzer, easy operation
- Powerful 400w halogen heater for quick and accurate measurements
- Compact design, backlit LCD graphic display
- Windows direct function for easy data transfer
- 4 different Measurement modes
- internal memory for 100 data sets and 10 programs incl. sample description
- RS 232 and USB interface
- 3 years warranty



Capacity	60 g
Readability	0.001 g / 0,01%
Reproducibility	0.002 g / 0,05%
Pan size	95 mm
Measuring units	Gram, %Moisture, % solids
Stabilisation time	2-3 sec
Internal Memory	10 Programs, 100 Test results
Calibration	automatic external
Heating source	400 W Halogen heater, 50° - 160°C
Drying programs	Standard, Drying in steps, gentle and fast
Display	huge, backlit LCD-Display