

## LabJacks LJ-Series

Brand new line of LabJacks. Designed with a height-adjustable platform, for supporting and adjusting the height of general lab devices such as hot plates, water baths, etc.; as well as glassware; mainly made of aluminum alloy materials.

The lab jacks are provided with good corrosion resistance and resistance to many chemicals.

### Features at a glance

- Easy adjustment with an adjusting wheel
- High stability and torsional rigidity
- stable stand
- High load carrying capacity and resilience
- Fully autoclavable at 121°C for sterile application
- 2 years warranty



LJ-01



LJ-02



LJ-03



Aluminum alloy



Adjusting wheel

Model	LJ-01	LJ-02	LJ-03
Base Size (mm)	80 x 76	140 x 120	150 x 121
Vertical Lift (mm)	48 x 148	58 x 242	78 x 278
Static Load (kg)	3,6	10	20
Dynamic Load (kg)	2	5,2	8

