

New Digital Infrared Vortexer RS-VA 20 and RS-VA 30



RS-VA 20

- Adjustable speed ranges from gentle to vigorous mixing to suit various needs
- LCD-screen displays real-time speed and time for easy adjustment and monitoring
- Stable with pressure-mounted accessories to prevent slippage and splashing
- Brushless DC motor for high stability
- Maintenance-free operation
- Sturdy base ensures stability during operation for safety

2 years warranty

RS-VA 30

- Compatible for a wide range of adapter for tubes and microplates
- Stable operation and sturdy construction
- Real time LCD display, Timer function
- Hands-free operation by IR-supported activation of the device
- Brushless motor, adjustable
- Wide speed range up to 3,000 rpm

Technical Data	RS-VA 20	RS-VA 30
Shaking Movement	orbital	
Amplitude	4 mm	4.5 mm
Speed	100 - 3,000 rpm	
Speed/Time adjustment	knob	
Display	LCD	
Operating Mode	IR, continous	Touch, continous
Time range	1 s - 99 min 59 s	
Dimension (LxWxH)	150 x 140 x 140 mm	205 x 165 x 90 mm
Protection class	IP 21	
Weight	3.4 kg	3.9 kg
Power Supply	100 - 240 V • 50/60 Hz	



RS-VA 30

Accessories for RS-VA 20 and RS-VA 30 on request

