

15-Place Stirrer RSM-03-15K

Newly developed 15-place magnetic stirrer is designed for series tests under reproducible conditions with stable performance, good reproducibility, high accuracy and easy operation.

RSM-03-15K

- 15 stirring positions are divided into 3 independent groups and each group can be controlled respectively at different speeds and time settings
- LCD displays show the current speeds and time
- Maintenance-free, brushless DC motor
- Aluminum laser engraving work plate surface, anti-abrasive, skid and corrosion proof
- Buzzer function - when the time countdown ends, the buzzer will remind the user by sound



RSM-03-15K

Technical Data	RSM-03-15K
Voltage • Frequency	100 - 240 V • 50/60 Hz
Power (Watt)	54
Stirring points	15
Max. stirring volume (Water)	15 x 0,4 l
Max. length of stirring bar	40 mm
Motortype	Brushless DC-motor
Speed range / Accuracy	100 - 1,500 / 10 rpm
Display	LCD
Plate dimension	460 x 304 mm
Timer function	1min - 99h 59min
Weight	9.2 kg
Ambient Temperature	5 - 40 °C
Relative Humidity	80 %
Safety class acc. DIN 60529	IP 21
Dimension (W x D x H)	615 x 307 x 90 mm